# PENIMASTER<sup>®</sup> PRO

de

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# Manual

en

PeniMaster® is a registered trademark of  
MSP Concept GmbH & Co. KG

Patent No. DE 10 2007 026 063

Utility model No. DE 20 2008 014 196

Utility model No. DE 20 2009 001 971

Utility model No. DE 20 2007 008 831

Utility model No. DE 20 2007 016 874

Utility model No. DE 20 2007 016 875

© Copyright 2018 MSP Concept GmbH & Co. KG

Last update: 02/2018. All rights reserved.

## Copyright

The present documentation is protected by copyright. All rights, particularly the right of reproduction and distribution, as well as translation, shall remain reserved.

No part of this publication must be reproduced in any way by photocopy, microfilm or other procedures, or be transferred into a language processible by machines, particularly data processing machines, without the written consent of the holder of the rights.

All rights of reproduction shall remain reserved.

Errors and modifications to the contents, as well as modifications serving technological progress, shall remain reserved.

# Parts list and designation of parts PeniMaster® PRO (Weight pulling force generator)

Art. No. ZPMP-007

en

- 1 PeniMaster® PRO Basic system**  
Art. No. PMP-001  
(not included)

- 2 Weight pulling force generator for PeniMaster® PRO system**  
Art. No. ZPMP-007

Base weight (250 g) with safety system (2.1) and screw-mountable stainless steel round weights (2.6). Pulling force adjustable in steps of 125 g (between 250 g and 1,250 g). Freely expandable with additional weights.

- 2.1 Base weight 250 g**  
Art. No. ZPMP-008  
Round weight with internal thread for screw-mounting of additional weights (2.6) and integral safety belt (2.3 – 2.5).
- 2.2 Coupling socket**  
Connects the weight pulling force generator (2) to PeniMaster® PRO (1).

**2.3 Safety cord**  
Tear-resistant, with fixed connection to the base weight.

**2.4 Shoulder belt**  
Wide, eudermic tissue belt. Worn like a sash around the shoulders.

**2.5 Ladder lock buckles**  
To adjust the length of the shoulder belt [2.4].

**2.6 Round weights**  
Round weights made of stainless steel. Screw-mounted to the base weight (250 g) [2.1] and to each other.

Included:  
2 x 125 g (ZPMP-009),  
1 x 250 g (ZPMP-010),  
1 x 500 g (ZPMP-011).

**3 Coupling retainer**  
Art. No. ZPMP-012  
Secures the connection of PeniMaster®PRO [1] and weight pulling force generator [2] by retaining the coupling socket [2.2].

## Application weight expander

### Before you start

The weight expander is an accessory for the PeniMaster®PRO basic system. Please therefore also read the PeniMaster®PRO basic system instructions on correct use and the risks and side effects.

### en Safety notice

The weights used always pose a kinetic risk (drop energy)!


- For safety reasons, always wear the safety belt and use the coupling retainer according to these instructions.
- Do not let the weight expander attached to the penis drop! There is a risk of injury if dropped just a few centimetres!
- Make slow and controlled movements when the weight expander is attached to the penis. Do not bounce or jump up and down. Walk slowly up and down stairs. Do not engage in any sporting activities.

### Preparation for application:

[!] PeniMaster®PRO already has to be applied to the penis (see page 20 et seq. of the instructions for PeniMaster®PRO).


### 1. Application as weight expander

- a) Slip one arm through the belt and pull it over your head. Worn on the body, the ladder lock buckle (2.5) must point towards the front.


a) Apply

- b) Using the ladder lock buckle (2.5), adjust the shoulder belt (2.4) so that its end (with the metal ring) is approximately at navel level.


*b) Adjust length using the ladder lock buckle*

- c) Press the narrow segment of the coupling head of PeniMaster®PRO (1) into the protrusion of the coupling socket (2.2) and pull on the base weight (2.1). In doing so, both parts are coupled together.

*c) Couple together*

- d) Open the coupling retainer (3) by pressing the wings together and clip the coupling retainer (3) on PeniMaster®PRO (1) as shown.


d) Clip on coupling retainer

- e) Increase the pulling force as needed by screw-mounting the appropriate round weights (2.6) to the base weight (2.1). (Note: You can also do this prior to step c) if you already know the needed pulling force).

The following pulling forces can be adjusted with the weight set (2.6):


- 250 g = base weight
- 375 g = base weight + 125 g
- 500 g = base weight + 250 g
- 625 g = base weight + 125 g + 250 g
- 750 g = base weight + 500 g
- 875 g = base weight + 125 g + 500 g
- 1000 g = base weight + 250 g + 500 g
- 1125 g = base weight + 125 g + 250 g + 500 g
- 1250 g = base weight + 125 g + 125 g + 250 g + 500 g

Note: The pulling force can be further increased by additional weights (not included).



e) Adjust pulling forces

Correctly applied weight expander system.

Variants

*Correctly applied*

# Legal notice

**Manufacturer of the penis expander system PeniMaster®PRO and publisher of this manual and these instructions for use**

MSP Concept GmbH & Co. KG  
Emser Str. 39d  
D-10719 Berlin  
Germany

Telephone: +49 (0)30 / 23 63 85 65  
Fax: +49 (0)30 / 21 75 26 11  
E-Mail: info@penimaster.de

District Court of Charlottenburg  
Commercial Register A No. 33128  
VAT ID No.: DE202064651

represented by

MSP Concept Verwaltungs GmbH  
Place of business: Berlin  
District Court of Charlottenburg  
Commercial Register B No. 83666

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**tr** MSP Concept Verwaltungs GmbH
ab-design GmbH, Berlin